

## Electroporation of *E. coli* cells

- Thaw cells on ice (~5 min)
- For each transformation:
  - Chill one 0.5 ml tube and a clean electroporation cuvette on ice.
  - Add 30  $\mu$ l of nanopure water and an appropriate amount (\*) of DNA to each tube.
  - Add 25  $\mu$ l of competent cells to each tube, mix by pipetting.
  - Transfer content to cuvette.
  - Before pulsing cells, tap cuvette to bring all liquid down, wipe metal sides with Kimwipe to remove water drops. Pulse with 1800 Volts. Time constant should be superior to 3.
  - Add 500  $\mu$ l of LB to each cuvette, mix, aspirate out and transfer to a 14 ml Falcon tube. Let recover for 30 min at 37°C with shaking.
  - Plate on LB agar plate supplemented with appropriate antibiotic. Incubate overnight at 37°C.

Testing competent cells.

Use pUC18 as a standard. Keep a stock at 10 ng/ $\mu$ l, dilute down by one thousand-fold to 10 pg/ $\mu$ l before transformation and use 1  $\mu$ l (10 picograms). Plate only one tenth of the final 556  $\mu$ l. Count the next day.

(\*) For regular plasmid transformation, one nanogram of DNA should give rise to about 100,000 transformants assuming an efficiency of  $10^8$  transformants per  $\mu$ g.

Note about ligations. Ligations contain salt. Salt will induce your electroporation to arc which will obliterate all your cells. At max, use 0.5  $\mu$ l of your ligation mixture per transformation. If more is required, perform an ethanol precipitation followed by a rinse with 70% ethanol prior to transformation.

### ***Washing cuvettes after use.***

Cuvettes can almost be reused indefinitely provided they're well washed.

After use, flush hard with distilled water in the sink for at least 30 secs per cuvette.

Add 0.2 M HCl to each cuvette, let stand for 10 min.

Flush hard with distilled water again.

Add 100% ethanol, let stand for 10 min. Empty, dry out inverted on your bench, store.